



Pevco Sentry™

- > Pevco Sentry™ provides healthcare facility engineers with automated email notification of pneumatic tube system transaction details, advisories and alarms.



As easy as email

Pevco Sentry™ is as easy as email so there is no need to train staff on new technology to reap its benefits. Pevco Sentry™ communicates with any device that is email-enabled – from desktop computer to mobile phone. Most major cellular service providers convert emails to text messages, so even a no-frills mobile phone can receive Sentry messages.

No additional technology investment required

Pevco Sentry™ is included standard with Pevco computer control centers equipped with Smart System software v2.28.0 or higher. If your system is running an earlier version of software, call Pevco's Help Desk at 800.296.7382 to schedule an upgrade. A Pevco service technician can be dispatched to your facility to put Sentry to work for you.

Implement system-wide or customize to fit your needs

Pevco Sentry™ can be customized to communicate with all levels of pneumatic tube system engineers and managers. By filtering emails to each addressee by station, and type of event, communication is streamlined. System engineers receive only the critical information they need to keep the system – and their facility – running smoothly.

Improve communication, productivity and response time with Pevco Sentry™

Here's how Sentry can work for you:

- Know when a station's bin is full and needs carriers removed to minimize downtime.
- Be alerted if a carrier fails to liftoff or fails to reach its destination to help avoid patient care delays or misplaced products.
- Receive notice of system-generated mechanical advisories and alarms so you can react quickly and reduce downtime, no matter where you are.

For more information about Pevco Sentry™, or other reliable Pevco products, go to www.pevco.com or call 800.296.7382.

